










# iOS HS SDK Document

## 1. Relevant file and Frameworks Configuration

### 1) Need to import HS SDK including:

HSHeader.h、HSMacroFile.h、HS3.h、HS3Controller.h、HS4.h、HS4Controller.h、HS5.h、HS5Controller.h、iHealthLibrary(x.x.x).a (Supports iOS 6.0 and above)

### 2) Frameworks

▼ Link Binary With Libraries (9 items)		×
Name	Status	
 MediaPlayer.framework	Required ⬆⬇⬆	
 SystemConfiguration.framework	Required ⬆⬇⬆	
 Accelerate.framework	Required ⬆⬇⬆	
 ExternalAccessory.framework	Required ⬆⬇⬆	
 CoreBluetooth.framework	Required ⬆⬇⬆	
 CoreGraphics.framework	Required ⬆⬇⬆	
 iHealthLibrary(1.0.7).a	Required ⬆⬇⬆	
 UIKit.framework	Required ⬆⬇⬆	
 Foundation.framework	Required ⬆⬇⬆	
+ — Drag to reorder frameworks		

### 3) Configuration

Add an “item” in “Info”

Add one new item in “Supported External accessory protocols”: com.ihealth.sc221

Add two new items in “Required background modes”: App communicates with an accessory, App communicates using Core Bluetooth

▼ Custom iOS Target Properties		
Key	Type	Value
Bundle versions string, short	String	1.0
Bundle identifier	String	com.zhang.xxx.\${PRODUCT_NAME:rfc1034identifie
InfoDictionary version	String	6.0
Main storyboard file base name	String	Main
Bundle version	String	1.0
▼ Required background modes	Array	(2 items)
Item 0	String	App communicates using CoreBluetooth ⬆⬇⬆
Item 1	String	App communicates with an accessory ⬆⬇⬆
▼ Supported external accessory protocols	Array	(6 items)
Item 0	String	com.jiuan.P930
Item 1	String	com.jiuan.BPV21
Item 2	String	com.jiuan.BPV20
Item 3	String	com.ihealth.sc221
Item 4	String	com.jiuan.BGV30
Item 5	String	com.jiuan.BGV31

## 2. Operation Procedure

### 1) HS3 instructions

Register plug-in device HS3 info: HS3ConnectNoti

Initialize HS3 controller class:

```
HS3Controller*controller=[HS3Controller sharedInstance];
```

Access control class instance after receiving HS3Controller:

```
NSArray *hsDeviceArray = [controller getAllCurrentHS3Instance];
```

```
HS3 *hsInstance = [hsDeviceArray objectAtIndex:i];
```

Use hsInstance to call HS3 related communication methods.

## 2) HS4 process instruction

Register plug-in device HS4 info: HS4ConnectNoti

Initialize HS4 controller class:

```
HS4Controller*controller=[HS4Controller sharedInstance];
```

Access control class instance after receiving HS4Controller:

```
NSArray *hsDeviceArray = [controller getAllCurrentHS4Instance];
```

```
HS4 *hsInstance = [hsDeviceArray objectAtIndex:i];
```

Use hsInstance to call HS4 related communication methods.

## 3) HS5 process instruction

Register plug-in device HS5 info: HS5ConnectNoti

Initialize HS5 controller class:

```
HS5Controller*controller=[HS5Controller sharedInstance];
```

Access control class instance after receiving HS5Controller:

```
NSArray *hsDeviceArray = [controller getAllCurrentHS5Instance];
```

```
HS5 *hsInstance = [hsDeviceArray objectAtIndex:i];
```

Use hsInstance to call HS5 related communication methods.

# 3. HS3 Interface Instructions

## 1) Establish memory and measurement connection

-(void)commandInitWithUser:(User \*)tempUser

Authentication:(BlockUserAuthentication)disposeAuthenticationBlock

TransferMemoryData:(StartHS3Transmission)startTransmission

UploadDataNum:(DisposeHS3UploadDataNum)hs3UploadDataNum

DisposeProgress:(DisposeHS3Progress)progress

MemoryData:(HS3MemoryData)memoryData

FinishTransmission:(FinishHS3Transmission)finishTransmission

StableWeight:(StableHS3Weight)stableWeight

DisposeErrorBlock:(DisposeHS3ErrorBlock)disposeHS3ErrorBlock;

### Import Parameters:

tempUser included properties: clientID, clientSecret, userID, and height

userID: either email or mobile phone number (mobile phone number is not currently supported yet)

height: the height of a user (cm);

clientID & clientSecret: the only identification for users of the SDK, requires registration from iHealth administrator, please email: [lvjincan@ihealthlabs.com.cn](mailto:lvjincan@ihealthlabs.com.cn) for more information.

**Return Parameters:**

disposeAuthenticationBlock: The return parameters of 'userid', 'height', 'clientID', and 'clientSecret' after verification

**The interpretation for the verification:**

UserAuthen\_RegisterSuccess: New-user registration succeeded.

UserAuthen\_LoginSuccess: User login succeeded.

UserAuthen\_CombinedSuccess: The user is an iHealth user as well, measurement via SDK has been activated, and the data from the measurement belongs to the user.

UserAuthen\_TrySuccess: Testing without internet connection succeeded.

UserAuthen\_InvalidateUserInfo: Userid/clientID/clientSecret verification failed.

UserAuthen\_SDKInvalidateRight: SDK has not been authorized.

UserAuthen\_UserInvalidateRight: User has not been authorized.

UserAuthen\_InternetError: Internet error, verification failed.

The measurement via SDK will be operated in the case of 1-4, and will be terminated if any of 5-8 occurs. The interface needs to be re-called after analyzing the return parameters. Notice: when a new user registers via SDK, an 'iHealth disclaimer' will pop up automatically, and will require the user to agree in order to continue. SDK applications require an Internet connection; there is 10-day trial period if the SDK cannot connect to the internet, the SDK is fully functional during tryout period, but will be terminated without a working internet connection after 10 days.

startTransmission: Start Memory transmission, Success: Yes, Fail: No.

hs3UploadDataNum: Memory Number, 0~200.

progress: Memory transmission progress, 0.0~1.0.

memoryData: Record data including weight(kg), measure time, coordinated key: weight, date.

finishTransmission: Finish memory transmission.

stableWeight: Stable weight (Kg)

disposeHS3ErrorBlock: Error code in measurement process

error code definition: refer to 'error' in Section 6: HS3 error instructions.

## **(2) Turn off Bluetooth Connection**

This method can be called only for hsInstance.HS3 with FirmwareVersion>=1.0.2

-(void)commandTurnOffBTConnectAutoResult:(DisposeResult)disposeResult

DisposeErrorBlock:(DisposeHS3ErrorBlock)disposeHS3ErrorBlock;

**Return parameters:**

disposeResult: YES' means measurement has been terminated, 'NO' means termination failed.

disposeHS3ErrorBlock: refer "error" in Section6: HS3 error instruction.

### (3) Turn on Bluetooth Connection

This method can be called only for hslInstance.HS3 with FirmwareVersion $\geq$ 1.0.2

-(void)commandTurnOnBTConnectAutoResult:(DisposeResult)disposeResult

DisposeErrorBlock:(DisposeHS3ErrorBlock)disposeHS3ErrorBlock;

#### **Return parameters:**

disposeResult: 'YES' means measurement has been terminated, 'NO' means termination failed.

disposeHS3ErrorBlock: refer to "error" in Section 6: HS3 error instruction.

## 4. HS4 Interface Method Instruction

### 1) Establish memory and measurement connection

When using the SDK for the first time, measuring method needs to be called to finish user verification.

-(void)commandMeasureWithUnit:(HSUnit)tempUnit andUser:(User \*)tempUser

Authentication:(BlockUserAuthentication)disposeAuthenticationBlock

Weight:(UnStableWeight)unStableWeight StableWeight:(StableWeight)stableWeight

DisposeErrorBlock:(DisposeHS4ErrorBlock)disposeErrorBlock;

#### **Import parameter:**

tempUnit: Unit displayed on HS4: HSUnit\_Kg、HSUnit\_LB、HSUnit\_ST。

tempUser: Properties included: clientID, clientSecret, userID, height.

userID、clientID、clientSecret, see the instruction of HS3

#### **Return parameters:**

disposeAuthenticationBlock, see the instruction of HS3

unStableWeight: Current weight, (Kg)

stableWeight: Stable weight, (Kg)

disposeErrorBlock: error code

#### **Error code definition:**

refer to "error" in Section 6: HS4 error instruction.

### 2) Upload memory data

If user doesn't pass the verification, HS4UserInvalidate will be returned for calling this method, user information invalid.

-(void)commandTransferMemoryWithUser:(User \*)tempUser

memoryData:(StartHS4Transmission)startTransmission

DisposeProgress:(DisposeProgress)progress

MemoryData:(MemoryData)memoryData

FinishTransmission:(FinishHS4Transmission)finishTransmission

DisposeErrorBlock:(DisposeHS4ErrorBlock)disposeErrorBlock;

#### **Import parameters:**

tempUser, included properties: userID, refer to the instructions for HS3

#### **Return parameters:**

startTransmission: Start Memory transmission.

progress: Memory transmission progress, 0.0~1.0.

memoryData:Record data including weight (kg), measurement time, coordinated key: weight, date.

finishTransmission: Finish memory transmission.

disposeErrorBlock: Record the error code in uploading process.

Error code definition: refer to "error" in Section 6: HS4 error instruction.

## 5. HS5 Interface Method Instruction

### 1) Establish memory and measurement connection

-(void)commandCreateUserManageConnectWithUser:(User \*)tempUser

Authentication:(BlockUserAuthentication)disposeAuthenticationBlock

currentUserSerialNub:(CurrentSerialNub)serialNub

deviceUserList:(MemoryUserListHS5Data) MemoryUserListHS5Data

Disposehs5ErrorBlock:(DisposeHS5ErrorBlock)disposeErrorBlock;

#### Import Parameters:

tempUser, included properties: userID、clientID、clientSecret。

userID、clientID、clientSecret、disposeAuthenticationBlock, refer to the instructions of HS3.

#### Return Parameters:

disposeAuthenticationBlock: **UserAuthen\_TrySuccess** is invalidate for HS5.

SerialNub: The only identification of a user, should be reserved in third party apps, provide to other Apps

MemoryUserListHS5Data: Existing user info in HS5, including serialNub、Position of users. Related key: serialNumber、position

disposeErrorBlock: error codes in transmission process, refer to Section 6: Errors in HS5

### 2) Create new user

Use the function if the SerialNub of current user is not included in user list of HS5 and the user number is less than 20

-(void)commandCreateUser:(User \*)tempUser position:(uint8\_t)tempPosition

DisposeHS5Result:(DisposeHS5Result)disposeHS5Result

Disposehs5ErrorBlock:(DisposeHS5ErrorBlock)disposeErrorBlock;

#### Import parameters:

tempUser, included properties: serialNub、birthday、height、isAthlete、sex。

serialNub: The only identification for user

birthday: User's birthday, NSDate

height: User's height, (cm)

isAthlete: If user is an athlete, UserIsAthlete\_No: no, UserIsAthlete\_Yes: Yes

sex: User's sex, UserSex\_Female: female, UserSex\_Male: male

tempPosition: Position of user, range: 0~19, pick from empty position that has never been occupied

#### Return Parameters:

disposeHS5Result: Success:Yes, Fail:No.

disposeErrorBlock: error codes in transmission process, refer to Section 6: Errors in HS5

### 3) Edit User Info

If the SerialNub of current user already exists in the user list of HS5, this will be used when user info is changed.

```
-(void)commandModifyUser:(User *)tempUser  
DisposeHS5Result:(DisposeHS5Result)disposeHS5Result  
Disposehs5ErrorBlock:(DisposeHS5ErrorBlock)disposeErrorBlock;
```

#### Import Parameters:

tempUser, included properties: serialNub、birthday、height、isAthlete、sex, instructions refer to “new user Api”.

#### Return Parameters:

disposeHS5Result: Success:Yes, Fail:No.  
disposeErrorBlock: error codes in transmission process, refer to Section 6: Errors in HS5

### 4) Delete Specified User

This will be used deleting account in HS5.

```
-(void)commandDelteUser:(User *)tempUser  
DisposeHS5Result:(DisposeHS5Result)disposeHS5Result  
Disposehs5ErrorBlock:(DisposeHS5ErrorBlock)disposeErrorBlock;
```

#### Import Parameters:

tempUser, included properties: serialNub, instructions refer to “new user Api”.

#### Return Parameters:

disposeHS5Result: Success:Yes, Fail:No.  
disposeErrorBlock: error codes in transmission process, refer to Section6: Errors in HS5

### 5) Establish Measurement Connection

```
-(void)commandCreateMeasureWithUser:(User *)tempUser  
unStableWeight:(UnStableHS5Weight)unStableHS5Weight  
MeasureWeight:(StableHS5Weight)stableHS5Weight  
ImpedanceType:(ImpedanceWeight)impedanceWeight  
BodyCompositionMeasurements:(BodyCompositionMeasurements)bodyCompositionMe  
asurements Disposehs5ErrorBlock:(DisposeHS5ErrorBlock)disposeErrorBlock;
```

#### Import Parameters:

tempUser, included properties: serialNub、height, instructions refer to “new user Api”.

#### Return Parameters:

unStableHS5Weight: Current weight, (kg)  
stableHS5Weight: Stable weight, (kg)  
impedanceWeight: Weight by impedance, (kg)  
bodyCompositionMeasurements: body info, includes weight(kg), fat content(%), water content(%), muscle content(%), bone mass, visceral fat level, DCI(Kcal). keys: weight, weightFatValue, waterValue, muscleValue, skeletonValue, VFatLevelValue, DCIValue  
disposeHS5Result: Success:Yes, Fail:No.  
disposeErrorBlock: error codes in transmission process, refer to Section 6: Errors in HS5

## 6) create memory upload connection

```
-(void)commandCreateMemoryWithUser:(User *)tempUser uploadConnect:  
(DisposeHS5Result)disposeHS5Result  
TransferMemoryData:(StartHS5Transmission)startTransmission  
DisposeProgress:(DisposeHS5Progress)progress  
MemoryData:(MemoryHS5Data)memoryData  
FinishTransmission:(FinishHS5Transmission)finishTransmission  
Disposehs5ErrorBlock:(DisposeHS5ErrorBlock)disposeErrorBlock;
```

### Import Parameters

tempUser, included properties: serialNub, instructions refer to “new user Api”.

### Return Parameters:

disposeHS5Result: Success:Yes, Fail:No.

startTransmission: Start memory transmission.

progress: Memory transmission progress, 0.0~1.0.

memoryData: Record data, More details and key refer Measure API. Additionally add time-measure property, related key: date.

finishTransmission: Finish Memory Transmission

disposeErrorBlock: error codes in transmission process, refer to Section 6: Errors in HS5

## 7) Acquire HS5 data from cloud

When connected to the internet, HS5 could upload measurement to cloud to be recorded.

Get data in cloud using the below API.

```
[[HS5Controller sharedInstance] getDownloadDataFromCloud:];  
-(void)getDownloadDataFromCloud:(BlockDataFromCloud )dataBlock;
```

### Return Parameters:

dataBlock: the collection of measurement from cloud, related key refer to Measure API.

## 5. Related supplementary parameter instruction

Device Connection info: **HS3ConnectNoti**、**HS4ConnectNoti**、**HS5ConnectNoti**

Device disconnection info: **HS3DisConnectNoti**、**HS4DisConnectNoti**、

**HS5DisConnectNoti**

When connecting to multiple devices, separate devices by deviceId

## 6. Error code instruction

(1)HS3:

```
typedef enum{  
    HS3DeviceLowPower = 1, //Low battery  
    HS3DeviceEr2, //Weight capacity is exceeded  
    HS3DeviceEr4, //The Scale calibration error
```

```

HS3DeviceEr7, //Movement while measuring
HS3DataZeor, //No memory
HS3DeviceDisconnect, //Device disconnect
HS3DeviceSendTimeout, //Communication error
HS3UserInvalidate = 111//User verify error
}HS3DeviceError;

```

(2)HS4:

```

typedef enum{
HS4DeviceLowPower = 1, // Battery level is low
HS4DeviceEr0, // The Scale failed to initialize
HS4DeviceEr1, // Maximum weight has been exceeded
HS4DeviceEr2, // The Scale can't capture a steady reading
HS4DeviceEr4, // Bluetooth connection error
HS4DeviceEr7, // Movement while measuring
HS4DeviceEr8, //Invalidate
HS4DeviceEr9, // Scale memory access error
HS4DataZeor, // No memory
HS4DeviceDisconnect, //Device disconnect
HS4DeviceSendTimeout, // Communication error
HS4DeviceRecWeightError, //
HS4UserInvalidate = 111//User verify error
}HS4DeviceError;

```

(3)HS5:

```

typedef enum{
IHSCOverTimeError = 0, // Communication error
IHSCUserInScale=5, // Communication Error
IHSCLowPower=6, // Make sure batteries are installed correctly, if the problem persists,
replace with a new set of batteries.
IHSCScaleEr0=7, // Remove the batteries, wait 1 minute and then replace with a new set
of batteries.
IHSCScaleEr1=8, // The current weight may be beyond the measurement range of 330
lbs/150 kg.
IHSCScaleEr2=9, // Stand still on all four electrodes with bare feet.
IHSCScaleEr7=10, // Communication Error
IHSCScaleEr8=11, // Communication Error
IHSCScaleEr9=12, // Communication Error
IHScaleBusy=13, // Scale is busy
HS5DataZeor=14, // No memory
IHSCScaleCreateUserTestError,
HS5Disconnect, //Device disconnect
HS5UserInvalidate = 111//User verify error
}HS5DeviceError;

```



## 7. Demo