










Instruction Manual of BP SDK

1. Relevant file and Frameworks configuration

- A. **BPSDK, including:** BPHeader.h, BPMacroFile.h, BP3.h, BP3Controller.h, BP5.h, BP5Controller.h, BP7.h, BP7Controller.h, ABI.h, ABIController.h, BP3L.h, BP3LController.h, BP7S.h, BP7SController.h, KN550BT.h, KN550BTController.h, KD926.h, KD926Controller.h, iHealthLibrary(x.x.x).a, support iOS6.0+

B. Frameworks

▼ Link Binary With Libraries (9 items) x

Name	Status
 MediaPlayer.framework	Required ⇅
 SystemConfiguration.framework	Required ⇅
 Accelerate.framework	Required ⇅
 ExternalAccessory.framework	Required ⇅
 CoreBluetooth.framework	Required ⇅
 CoreGraphics.framework	Required ⇅
 iHealthLibrary(1.0.7).a	Required ⇅
 UIKit.framework	Required ⇅
 Foundation.framework	Required ⇅

+ — Drag to reorder frameworks

C. Configuration

Add new 'Item' in 'Info':

Add 2 new Item in '**Supported external accessory protocols**':
[com.jiuan.BPV20](#), [com.jiuan.P930](#), [com.jiuan.BPV21](#), [com.jiuan.BPV24](#),
[com.jiuan.BPAV10](#), [com.jiuan.BPV22](#),

Add 1 new Item in '**Required background modes**': App communicates with an accessory、App communicates using CoreBluetooth

▼ Custom iOS Target Properties

Key	Type	Value
Bundle versions string, short	String	1.0
Bundle identifier	String	com.zhang.xxx.\${PRODUCT_NAME:rfc1034identifie
InfoDictionary version	String	6.0
Main storyboard file base name	String	Main
Bundle version	String	1.0
▼ Required background modes	Array	(2 items)
Item 0	String	App communicates using CoreBluetooth
Item 1	String	App communicates with an accessory
▼ Supported external accessory protocols	Array	(6 items)
Item 0	String	com.jiuan.P930
Item 1	String	com.jiuan.BPV21
Item 2	String	com.jiuan.BPV20
Item 3	String	com.ihealth.sc221
Item 4	String	com.jiuan.BGV30
Item 5	String	com.jiuan.BGV31

2. Operation Procedure

(1) Operation procedure for BP3

- register plug-in device info: `BP3ConnectNoti;`
- initialize controller classes:

```
BP3Controller *controller = [BP3Controller
shareBP3Controller];
```

- Access control class instance after receive `BP3ConnectNoti:`

```
NSArray *bpDeviceArray = [controller
getAllCurrentBP3Instance];
BP3 *bpInstance = [bpDeviceArray objectAtIndex: 0];
```

- Using 'bpInstance' call communication module of the device

(2) Operation procedure for BP5

- register plug-in device info: `BP5ConnectNoti:`
- initialize controller classes:

```
BP5Controller *controller = [BP5Controller
shareBP5Controller];
```

- Access control class instance after receive `BP5ConnectNoti:`

```
NSArray *bpDeviceArray = [controller
getAllCurrentBP5Instance];
BP5 *bpInstance = [bpDeviceArray objectAtIndex: i];
```

- Using 'bpInstance' call communication module of the device

(3) Operation procedure for BP7

- register plug-in device info: `BP7ConnectNoti;`
- initialize controller classes:

```
BP7Controller *controller = [BP7Controller
shareBP7Controller];
```

- Access control class instance after receive `BP7ConnectNoti:`

```
NSArray *bpDeviceArray = [controller
getAllCurrentBP7Instance];
BP7 *bpInstance = [bpDeviceArray objectAtIndex: i];
```

- d) Using 'bpInstance' call communication module of the device

(4) Operation procedure for ABI

For ABI Measure (both arm and leg)

- a) Register plug-in device info: `ABISConnectNoti`;
b) Initialized controller class:
`ABIController *controller = [ABIController
shareABIController];`
c) Access controller class instance after receive `ABISConnectNoti`:
`ABI *bpInstance = [controller getCurrentABIInstace];`
d) Using 'bpInstance' call communication module of the device.

For Arm Measure (arm only)

- a) Register plug-in device info: `ArmConnectNoti`;
b) Initialized controller class:
`ABIController *controller = [ABIController
shareABIController];`
c) Access controller class instance after receive `ArmConnectNoti`:
`ABI *bpInstance = [controller getCurrentArmInstace];`
d) Using 'bpInstance' call communication module of the device.

(5) Operation procedure for BP3L

- a) register plug-in device info: `BP3LConnectNoti`;
b) initialize controller classes:
`BP3LController *controller = [BP3LController
shareBP3LController];`
c) Access control class instance after receive `BP3LConnectNoti`:
`NSArray *bpDeviceArray = [controller
getAllCurrentBP3LInstace];`
`BP3L *bpInstance = [bpDeviceArray objectAtIndex: 0];`
Using 'bpInstance' call communication module of the device

(6) Operation procedure for BP7S

- a) register plug-in device info: `BP7SConnectNoti`;
b) initialize controller classes:
`BP7SController *controller = [BP7SController
shareBP7SController];`
c) Access control class instance after receive `BP7SConnectNoti`:
`NSArray *bpDeviceArray = [controller
getAllCurrentBP7SInstace];`
`BP7S *bpInstance = [bpDeviceArray objectAtIndex: 0];`
Using 'bpInstance' call communication module of the device

(7) Operation procedure for KN550BT

a) register plug-in device info: `KN550BTConnectNoti`;

b) initialize controller classes:

```
KN550BTController *controller = [KN550BTController  
shareKN550BTController];
```

c) Access control class instance after receive `KN550BTConnectNoti`:

```
NSArray *bpDeviceArray = [controller  
getAllCurrentKN550BTInstance];  
KN550BT *bpInstance = [bpDeviceArray objectAtIndex: 0];
```

Using 'bpInstance' call communication module of the device

(8) Operation procedure for KD926

a) register plug-in device info: `KD926ConnectNoti`;

b) initialize controller classes:

```
KD926Controller *controller = [KD926Controller  
shareKD926Controller];
```

c) Access control class instance after receive `KD926ConnectNoti`:

```
NSArray *bpDeviceArray = [controller  
getAllCurrentKD926Instance];  
KD926 *bpInstance = [bpDeviceArray objectAtIndex: 0];
```

Using 'bp Instance' call communication module of the device

3. Interface Method:

BP3:

a) Establish measurement connection

```
-(void)commandStartMeasureWithUser:(NSString *)userID  
clientID:(NSString *)clientID clientSecret:(NSString  
*)clientSecret  
Authentication:(BlockUserAuthentication)disposeAuthentication  
Block pressure:(BlockPressure)pressure  
xiaoboWithHeart:(BlockXiaoaboWithHeart)xiaobo  
xiaoboNoHeart:(BlockXiaoaboNoHeart)xiaoboNoHeart  
result:(BlockMeasureResult)result errorBlock:(BlockError)error;
```

Import parameters:

'userID', the only identification for the user, by the form of email or cell phone #(cell-phone-# form is not supported temperately)

'clientID' and 'clientsecret' are the only identification for user of SDK, are required registration from iHealth administrator, please email: lvjincan@ihealthlabs.com.cn for more information.

Return parameter:

disposeAuthenticationBlock is the return parameter of 'userid', 'clientID', 'clientSecret' after the verification. The interpretation for the verification:

1. UserAuthen_RegisterSuccess, New-user registration succeeded.
2. UserAuthen_LoginSuccess, User login succeeded.
3. UserAuthen_CombinedSuccess, The user is iHealth user as well, measurement via SDK has been activated, and the data from the measurement belongs to the user.
4. UserAuthen_TrySuccess, Testing without Internet connection succeeded.
5. UserAuthen_InvalidateUserInfo, Userid/clientID/clientSecret verification failed.
6. UserAuthen_SDKInvalidateRight, SDK has not been authorized.
7. UserAuthen_UserInvalidateRight, User has not been authorized.
8. UserAuthen_InternetError, Internet error, verification failed.

The measurement via SDK will be operated in the case of 1-4, and will be terminated if any of 5-8 occurs. The interface needs to be re-called after analyzing the return parameters.

Notice: by the first time of new user register via SDK, 'iHealth disclaimer' will pop up automatically, and require the user agrees to continue. SDK application requires Internet connection; there is 10-day tryout if SDK cannot connect Internet, SDK is fully functional during tryout period, but will be terminated without verification through Internet after 10 days.

'Pressure': Pressure value in the process of measurement, the unit is 'mmHg'.

- **'Xiaobo'**: Wavelet data set including pulse rate
- **'XiaoboNoHeart'**: Wavelet data set without pulse rate
- **'Result'**: result of the measurement, including systolic pressure, diastolic pressure, pulse rate and irregular judgment. Relevant key: **time, sys, dia, heartRate, irregular**

· 'Error': refer 'error' in in part 5

b) Measurement termination

```
-(void)stopBPMeasure:(BlockStopResult)result  
ErrorBlock:(BlockError)error;
```

Return parameter:

- 'Result': 'YES' means measurement has been terminated, 'NO' means termination failed.
- 'error': refer 'error' in 'Establish measurement connection'

c) Battery life query

```
-(void)commandEnergy:(BlockEnergyValue)energyValue  
errorBlock:(BlockError)error;
```

Return parameter:

- energyValue: battery percentage, '80' stands for 80%.
- 'error': refer 'error' in 'Establish measurement connection'

BP5:

a) Establish measurement connection

```
-(void)commandStartMeasureWithUser:(NSString *)userID  
clientID:(NSString *)clientID clientSecret:(NSString  
*)clientSecret  
Authentication:(BlockUserAuthentication)disposeAuthentication  
Block pressure:(BlockPressure)pressure  
xiaoboWithHeart:(BlockXiaoaboWithHeart)xiaobo  
xiaoboNoHeart:(BlockXiaoaboNoHeart)xiaoboNoHeart  
result:(BlockMesureResult)result errorBlock:(BlockError)error;
```

Import parameters:

'userID', 'clientID', 'clientSecret',: refer in 'BP3'

Return parameter:

- 'disposeAuthenticationBlock': refer in 'BP3'
- 'Pressure': Pressure value in the process of measurement, the unit is 'mmHg'.
- 'Xiaoabo': Wavelet data set including pulse rate
- 'XiaoaboNoHeart': Wavelet data set without pulse rate
- 'Result': result of the measurement, including systolic pressure, diastolic pressure, pulse rate and irregular judgment. Relevant key: **time, sys, dia, heartRate, irregular**

· 'error': refer 'error' in in part 5

BP7:

a) Establish measurement connection

```
-(void)commandStartGetAngleWithUser:(NSString *)userID  
clientID:(NSString *)clientID clientSecret:(NSString  
*)clientSecret  
Authentication:(BlockUserAuthentication)disposeAuthentication  
Block angle:(BlockAngle)angleInfo  
errorBlock:(BlockError)error;
```

Import parameters:

'userID', 'clientID', 'clientSecret',: refer in 'BP3'

Return parameter:

'disposeAuthenticationBlock': refer in 'BP3'
'AngleInfo': angle info, including angle、isLeftHand.
'angle' stands for the inclination angle of the device, from 0 to 180; 'isLeftHand' stands for the judgment if the device is on the left hand, 'true' means left hand, otherwise is right hand.
'error': refer 'error' in in part 5

b) Measurement start:

When 'angle' is between 10-30, measurement starts:

Return parameter:

```
-(void)commandStartMeasure:(BlockPressure)pressure  
xiaoboWithHeart:(BlockXiaoaboWithHeart)xiaobo  
xiaoboNoHeart:(BlockXiaoaboNoHeart)xiaoboNoHeart  
result:(BlockMesureResult)result errorBlock:(BlockError)error;
```

- 'Pressure': Pressure value in the process of measurement, the unit is 'mmHg'.
- 'Xiaobo': Wavelet data set including pulse rate
- 'XiaoaboNoHeart': Wavelet data set without pulse rate
- 'Result': result of the measurement, including systolic pressure, diastolic pressure, pulse rate and irregular judgment. Relevant key: time, sys, dia, heartRate, irregular
- 'error': refer 'error' in in part 5

BP5 & BP7(self-adaptive)

a) Hypogenous query

-(void)commandFunction:(BlockDeviceFunction)founction
errorBlock:(BlockError)error;

Return parameter:

Function: judge if the device supports BT auto-connection, offline detection, and if the function on or off, corresponding KEY as haveBlue, haveOffline, blueOpen, offlineOpen. 'True' means yes or on, 'False' means no or off.

'error': refer 'error' in in part 5

b) set up BT auto-connection

-(void)commandSetBlueConnect:(BOOL)open
respond:(BlockBlueSet)blockBuleSet
errorBlock:(BlockError)error

Import parameter:

Open: True means on; False means off.

Return parameter:

blockBuleSet: device current statues, True means on, False means off.

'error': refer 'error' in in part 5

c) set up offline detection

-(void)commandSetOffline:(BOOL)open
respond:(BlockOfflineSet)blockOfflineSet
errorBlock:(BlockError)error

Import parameter:

Open: True means on; False means off.

Return parameter:

blockOfflineSet: device current statues, True means on, False means off.

'error': refer 'error' in in part 5

d) Battery life query

-(void)commandEnergy:(BlockEnergyValue)energyValue
errorBlock:(BlockError)error

Return parameter:

energyValue: battery percentage, '80' stands for 80%.

'error': refer 'error' in in part 5

e) Measurement termination

-(void)stopBPMeasureErrorBlock:(BlockStopSuccess)success
errorBlock:(BlockError)error

Return parameter:

'Result': 'YES' means measurement has been terminated, 'NO' means termination failed.

'error': refer 'error' in in part 5

f) upload offline data

-(void)commandBatchUpload:(BlockBatchCount)totalCount
progress:(BlockBatchProgress)progress
dataArray:(BlockBatchArray)uploadDataArray
errorBlock:(BlockError)error;

Return parameter:

TotalCount: item quantity of total data

Progress: upload completion ratio , from 0.0 to 1.0 or 0%~100%, 100% means upload completed

UploadDataArray:offline data set, including measurement time, systolic pressure, diastolic pressure, pulse rate, irregular judgment. corresponding KEY as time, sys, dia, heartRate, irregular

'error': refer 'error' in in part 5

ABI:

a) Establish measurement connection

-(void)commandStartMeasureWithUser:(NSString *)userID
clientId:(NSString *)clientId clientSecret:(NSString *)clientSecret
Authentication:(BlockUserAuthentication)disposeAuthentication
Block armPressure:(BlockPressure)armPressure
legPressure:(BlockPressure)legPressure
armXiaoBoWithHeart:(BlockXiaoBoWithHeart)armXiaoBo
legXiaoBoWithHeart:(BlockXiaoBoWithHeart)legXiaoBo
armXiaoBoNoHeart:(BlockXiaoBoNoHeart)armXiaoBoNoHeart
legXiaoBoNoHeart:(BlockXiaoBoNoHeart)legXiaoBoNoHeart
armResult:(BlockMeasureResult)armResult

legResult:(BlockMeasureResult)legResult
errorBlock:(BlockError)error;

Import parameters:

‘userID’, ‘clientID’, ‘clientSecret’,: refer in ‘BP3’

Return parameter:

‘disposeAuthenticationBlock’: refer in ‘BP3’
‘armPressure’:Upper-arm blood pressure value during measurement, unit as mmHg
‘legPressure’: Ankle blood pressure value during measurement, unit as mmHg.
‘armXiaobo’:Wavelet value of upper-arm BPM, with heartbeats.
‘legXiaobo’: Wavelet value of ankle BPM, with heartbeats.
‘armXiaoboNoHeart’: Wavelet value of upper-arm BPM, without heartbeats.
‘legXiaoboNoHeart’: Wavelet value of ankle BPM, without heartbeats.
‘armResult’: BP value of upper-arm BPM, including time, sys, dia, heartRate, irregular heartbeat.
‘legResult’: BP value of ankle BPM, including time, sys, dia, heartRate, irregular heartbeat.

‘error’: error codes details in part 5.

b) Stop measuring

-(void)stopABIMeasureErrorBlock:(BlockStopSuccess)success
errorBlock:(BlockError)error;

Return parameter:

‘Result’: ‘YES’ means measurement has been terminated, ‘NO’ means termination failed.
‘error’: refer ‘error’ in in part 5.

c) Check battery life

-(void)commandQueryEnergy:(BlockEnergyValue)armEnergy
leg:(BlockEnergyValue)legEnergy errorBlock:(BlockError)error;

Return parameters:

‘armEnergy’: Returns battery ratio of upper-arm BPM, 80 means

80%.

'legEnergy': Returns batter ratio of ankle BPM, 80 means 80%
'error': refer 'error' in in part 5.

d) Establish Arm measurement connection

```
-(void)commandStartMeasureWithUser:(NSString *)userID  
clientID:(NSString *)clientID clientSecret:(NSString  
*)clientSecret  
Authentication:(BlockUserAuthentication)disposeAuthentication  
Block armPressure:(BlockPressure)armPressure  
armXiaoboWithHeart:(BlockXiaoaboWithHeart)armXiaoabo  
armXiaoaboNoHeart:(BlockXiaoaboNoHeart)armXiaoaboNoHeart  
armResult:(BlockMesureResult)armResult  
errorBlock:(BlockError)error;
```

Import parameters:

See Api(1)

Return parameter:

See Api(a)

e) Stop Arm measuring

```
-(void)stopABIArmMeassureBlock:(BlockStopResult)result  
errorBlock:(BlockError)error;
```

Return parameter:

See Api(b)

f) Check Arm battery life

```
-(void)commandQueryEnergy:(BlockEnergyValue)armEnergy  
errorBlock:(BlockError)error;
```

Return parameters:

See Api(c)

BP3L:

a) Establish measurement connection

```
-(void)commandStartMeasureWithUser:(NSString *)userID  
clientID:(NSString *)clientID clientSecret:(NSString  
*)clientSecret  
Authentication:(BlockUserAuthentication)disposeAuthentication  
Block pressure:(BlockPressure)pressure  
xiaoboWithHeart:(BlockXiaoaboWithHeart)xiaobo  
xiaoboNoHeart:(BlockXiaoaboNoHeart)xiaoboNoHeart  
result:(BlockMesureResult)result errorCallback:(BlockError)error;
```

Import parameters:

'userID', the only identification for the user, by the form of email or cell phone #(cell-phone-# form is not supported temperately).

'clientID' and **'clientsecret'** are the only identification for user of SDK, are required registration from iHealth administrator, please email: lvjincan@ihealthlabs.com.cn for more information.

Return parameter:

'disposeAuthenticationBlock' is the return parameter of **'userid'**, **'clientID'**, **'clientSecret'** after the verification. The interpretation for the verification:

1. UserAuthen_RegisterSuccess,
New-user registration succeeded.
2. UserAuthen_LoginSuccess, User login succeeded.
3. UserAuthen_CombinedSuccess, The user is iHealth user as well, measurement via SDK has been activated, and the data from the measurement belongs to the user.
4. UserAuthen_TrySuccess, Testing without Internet connection succeeded.
5. UserAuthen_InvalidaterUserInfo,
Userid/clientID/clientSecret verification failed.
6. UserAuthen_SDKInvalidaterRight, SDK has not been authorized.
7. UserAuthen_UserInvalidaterRight, User has not been authorized.
8. UserAuthen_InternetError, Internet error, verification failed.

The measurement via SDK will be operated in the case of 1-4, and will be terminated if any of 5-8 occurs. The interface needs

to be re-called after analyzing the return parameters.

Notice: by the first time of new user register via SDK, 'iHealth disclaimer' will pop up automatically, and require the user agrees to continue. SDK application requires Internet connection; there is 10-day tryout if SDK cannot connect Internet, SDK is fully functional during tryout period, but will be terminated without verification through Internet after 10 days.

· **'Pressure'**: Pressure value in the process of measurement, the unit is 'mmHg'.

· **'Xiaobo'**: Wavelet data set including pulse rate

· **'XiaoboNoHeart'**: Wavelet data set without pulse rate

· **'Result'**: result of the measurement, including systolic pressure, diastolic pressure, pulse rate and irregular judgment. **Relevant key: sys, dia, heartRate, irregular**

· **'Error'**: refer 'error' in in part 5

b) Measurement termination

```
-(void)stopBPMeasure:(BlockStopResult)result
```

```
ErrorBlock:(BlockError)error;
```

Return parameter:

· **'Result'**: 'YES' means measurement has been terminated, 'NO' means termination failed.

· **'error'**: refer 'error' in 'Establish measurement connection'

c) Hypogenous query

```
-(void)commandFounction:(BlockDeviceFounction)founction
```

```
errorBlock:(BlockError)error;
```

Return parameter:

· **'Function'**: judge if the device supports BT auto-connection, offline detection, and if the function on or off, corresponding KEY as haveBlue, haveOffline, blueOpen, offlineOpen. 'True' means yes or on, 'False' means no or off.

· **'error'**: refer 'error' in in part 5.

d) Battery life query

```
-(void)commandEnergy:(BlockEnergyValue)energyValue
```

```
errorBlock:(BlockError)error;
```

Return parameter:

'energyValue': battery percentage, '80' stands for 80%.
'error': refer 'error' in in part 5.

BP7S:

a) Hypogenous query

```
-(void)commandFunction:(BlockDeviceFunction)founction  
errorBlock:(BlockError)error;
```

Return parameter:

- **'Function'**: judge if the device supports BT auto-connection, offline detection, and if the function on or off, corresponding KEY as haveAngleSensor、haveAngleSet、haveHSD、haveOffline、mutableUpload、selfUpdate. 'True' means yes or on, 'False' means no or off.
- **'error'**: refer 'error' in in part 5.

b) Battery life query

```
-(void)commandEnergy:(BlockEnergyValue)energyValue  
errorBlock:(BlockError)error;
```

Return parameter:

'energyValue': battery percentage, '80' stands for 80%.
'error': refer 'error' in in part 5.

c) upload offline data

```
-(void)commandTransferMemoryDataWithUser:(NSString  
)userID clientID:(NSString *)clientID  
clientSecret:(NSString *)clientSecret  
Authentication:(BlockUserAuthentication)disposeAuthenticat  
ionBlock totalCount:(BlockBachCount)totalCount  
pregress:(BlockBachProgress)progress  
dataArray:(BlockBachArray)uploadDataArray  
errorBlock:(BlockError)error;
```

Import parameters:

'userID', 'clientID', 'clientSecret',: refer in 'BP3'.

Return parameter:

‘totalCount’: item quantity of total data.
‘progress’:upload completion ratio , from 0.0 to 1.0 or 0%~100%, 100% means upload completed .
‘uploadDataArray’: offline data set, including measurement time, systolic pressure, diastolic pressure, pulse rate, irregular judgment. corresponding KEY as time、sys、dia、heartRate、irregular、hsdValue、chooseHand、startAngle、measureAngle. ‘error’: refer ‘error’ in in part 5.

KN550BT:

a) Hypogenous query

-(void)commandFunction:(BlockDeviceFunction)founction
errorBlock:(BlockError)error;

Return parameter:

- ‘Function’: judge if the device supports BT auto-connection, offline detection, and if the function on or off, corresponding KEY as haveBlue, haveOffline, blueOpen, offlineOpen. ‘True’ means yes or on, ‘False’ means no or off.
- ‘error’: refer ‘error’ in in part 5.

b) Battery life query

-(void)commandEnergy:(BlockEnergyValue)energyValue
errorBlock:(BlockError)error;

Return parameter:

‘energyValue’: battery percentage, ‘80’ stands for 80%.
‘error’: refer ‘error’ in in part 5.

c) upload offline data

-(void)commandTransferMemoryDataWithUser:(NSString *)userID clientID:(NSString *)clientID
clientSecret:(NSString *)clientSecret
Authentication:(BlockUserAuthentication)disposeAuthenticat
ionBlock totalCount:(BlockBachCount)totalCount
pregress:(BlockBachProgress)progress
dataArray:(BlockBachArray)uploadDataArray

errorBlock:(BlockError)error;

Import parameters:

‘userID’, ‘clientId’, ‘clientSecret’,: refer in ‘BP3’.

Return parameter:

‘totalCount’: item quantity of total data.

‘progress’: upload completion ratio, from 0.0 to 1.0 or 0%~100%, 100% means upload completed.

‘uploadDataArray’: offline data set, including measurement time, systolic pressure, diastolic pressure, pulse rate, irregular judgment. corresponding KEY as time、sys、dia、heartRate、irregular.

‘error’: refer ‘error’ in in part 5.

KD926:

a) Hypogenous query

-(void)commandFunction:(BlockDeviceFunction)function
errorBlock:(BlockError)error;

Return parameter:

- ‘Function’: judge if the device supports BT auto-connection, offline detection, and if the function on or off, corresponding KEY as haveBlue, haveOffline, blueOpen, offlineOpen. ‘True’ means yes or on, ‘False’ means no or off.

- ‘error’: refer ‘error’ in in part 5.

b) Battery life query

-(void)commandEnergy:(BlockEnergyValue)energyValue
errorBlock:(BlockError)error;

Return parameter:

‘energyValue’: battery percentage, ‘80’ stands for 80%.

‘error’: refer ‘error’ in in part 5.

c) upload offline data


```
-(void)commandTransferMemoryDataWithUser:(NSString *)userID clientID:(NSString *)clientID clientSecret:(NSString *)clientSecret Authentication:(BlockUserAuthentication)disposeAuthenticationBlock totalCount:(BlockBachCount)totalCount progress:(BlockBachProgress)progress dataArray:(BlockBachArray)uploadDataArray errorCallback:(BlockError)error;
```

Import parameters:

'userID', 'clientID', 'clientSecret',: refer in 'BP3'.

Return parameter:

'totalCount': item quantity of total data.

'progress' : upload completion ratio , from 0.0 to 1.0 or 0%~100%, 100% means upload completed .

'uploadDataArray': offline data set, including measurement time, systolic pressure, diastolic pressure, pulse rate, irregular judgment. corresponding KEY as time、sys、dia、heartRate、irregular.

'error': refer 'error' in in part 5.

4. Parameter supplementary instruction

Device connection status: BP3ConnectNoti, BP5ConnectNoti, BP7ConnectNoti, ABICConnectNoti, BP3LConnectNoti, BP7SConnectNoti, KN550BTConnectNoti, KD926ConnectNoti

Device disconnection status: BP3DisConnectNoti, BP5DisConnectNoti, BP7DisConnectNoti, ABIDisConnectNoti, ArmConnectNoti, ArmDisConnectNoti, BP3LDisConnectNoti, BP7SDisConnectNoti, KN550BTDisConnectNoti, KD926DisConnectNoti

'serialNumber' is for separating different device when multiple device have been connected.

5. Error

```
typedef enum {  
    BPErr0 = 0, //Unable to take measurements due to arm/wrist movements.  
    BPErr1, //Failed to detect systolic pressure  
    BPErr2, //Failed to detect diastolic pressure
```

iHealthLabs EUROPE

during inflation	BPErr3, //Pneumatic system blocked or cuff is too tight
during inflation	BPErr4, //Pneumatic system leakage or cuff is too loose
160 seconds	BPErr5, //Cuff pressure reached over 300mmHg BPErr6, //Cuff pressure reached over 15 mmHg for more than
	BPErr7, //Data retrieving error BPErr8, //Data retrieving error BPErr9, //Data retrieving error BPErr10, //Data retrieving error BPErr11, //Communication Error BPErr12, //Communication Error BPErr13, //Low battery BPErr14, //
199mmHg	BPErr15, //Systolic exceeds 260mmHg or diastolic exceeds
	BPErr16, //Systolic below 60mmHg or diastolic below 40mmHg BPErr17, //Arm/wrist movement beyond range BPNormalError = 30, //device error, error message displayed
automatically	BPOverTimeError, //Abnormal communication BPNoRespondError, //Abnormal communication BPBeyondRangeError, //device is out of communication range. BPDidDisconnect, //device is disconnected. BPAskToStopMeasure //measurement has been stopped.
	}BPDeviceError;

6. Demo