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▼ Link Binary With Libraries (9 items)		
Name	Status	
 MediaPlayer.framework	Required	⌵
 SystemConfiguration.framework	Required	⌵
 Accelerate.framework	Required	⌵
 ExternalAccessory.framework	Required	⌵
 CoreBluetooth.framework	Required	⌵
 CoreGraphics.framework	Required	⌵
 iHealthLibrary(1.0.7).a	Required	⌵
 UIKit.framework	Required	⌵
 Foundation.framework	Required	⌵
+ — Drag to reorder frameworks		

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```
-(void)commandStartP03MeasureData:(StartP03MeasureData)startP03MeasureData
Measure:(DisposeP03MeasureData)disposeP03MeasureData
FinishP03MeasureData:(FinishP03MeasureData)finishP03MeasureData
```

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finishP03MeasureData: 成功返回 YES, 失败返回 NO.

disposeErrorBlock: 见第 5 部分错误说明。

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-(void)commandDisposeP03DataCount:(DisposeP03DataCount)disposeP03DataCount

TransferMemoryData:(StartP03Transmission)startTransmission

Memory:(DisposeP03HistoryData)disposeP03HistoryData

DisposeP03WaveHistoryData:(DisposeP03WaveHistoryData)disposeP03WaveHistoryData

DisposeProgress:(DisposeP03ProgressData)progress

FinishTransmission:(FinishP03Transmission)finishTransmission

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-(void)commandQueryBatteryInfo:(DisposeP03Block)disposeBlock

DisposeErrorBlock:(DisposeP03ErrorBlock)disposeErrorBlock

DisposeBattery:(DisposeP03Battery)disposeBattery;

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```
typedef enum{  
    P03CommError = 0,      // Bluetooth Communication Error  
    P03AccessError,        // Flash (Data) Access Error  
    P03HardwareError,      // Irregular Hardware Error  
    P03PRbpmtestError,     // The SpO2 or pulse rate test result is beyond  
the measurement range of the system  
    P03UnknownError,       // Unknown Interference Detected  
    P03SendCommandFaild,   // Send failed  
    P03DeviceDisConect,    // Device is disconnected  
    P03DataZero,           // No data  
    P03UserInvalidate      // User authentication fails  
}P03ErrorID;
```

6. 测试 demo