

# iOS HS6 SDK Document







## 1. Relevant file and Frameworks Configuration

### 1) Need to import HS6 SDK including:

HS6User.h、iHealthHS6.h、HS6Library (x.x.x).a (Supports iOS 6.0 and above)

### 2) Frameworks

▼ Linked Frameworks and Libraries

Name	Status
 ExternalAccessory.framework	Required ⇅
 Accelerate.framework	Required ⇅
 HS6Library(1.0.2).a	Required ⇅
 SystemConfiguration.framework	Required ⇅
 Foundation.framework	Required ⇅
 UIKit.framework	Required ⇅
+ -	

## 2. Operation Procedure

### 1) Set Wifi for HS6

If set Wifi success, the Wifi symbol will become stable.

### 2) Link User and HS6

If link success, the data of HS6 will upload to cloud automatically, then you can download the data of user from cloud.

### 3) Unlink User and HS6 if needed

If unlink success, the data of HS6 will not upload to cloud.

## 3. Interface Instructions

### 1) Set HS6 Wifi

```
-(void)commandSetHS6WithPassWord:(NSString*)password  
disposeHS6SuccessBlock:(DisposeHS6SuccessBlock)disposeHS6SuccessBlock  
disposeHS6FailBlock:(DisposeHS6FailBlock)disposeHS6FailBlock  
disposeHS6EndBlock:(DisposeHS6EndBlock)disposeHS6EndBlock  
disposeHS6ErrorBlock:(DisposeHS6ErrorBlock)disposeHS6ErrorBlock;
```

#### Import Parameters:

password: Wifi password.

#### Return Parameters:

disposeHS6SuccessBlock: link success

disposeHS6FailBlock: user link failed

disposeHS6EndBlock: end set

disposeHS6ErrorBlock: error code

Error code definition:

refer to "error" : HS6 error instruction.

## (2) Link HS6

-(void)cloudCommandUserBinedQRDeviceWithUser:(HS6User \*)tempUser

deviceId:(NSString \*)tempDeviceID

BlockHS6UserAuthentication:(BlockHS6UserAuthentication)BlockHS6UserAuthentication

binedResult:(BinedQRDeviceBlock)result binedError:(BinedQRDeviceErrorBlock)error;

### Import parameter:

tempUser: Properties included: clientID, clientSecret, userID, height.

deviceId: device Mac. Scan the QR code to get the full Mac ID.

clientID & clientSecret: the only identification for users of the SDK, requires registration from iHealth administrator, please email: lvjincan@ihealthlabs.com.cn.com for more information.

### Return parameters:

BlockHS6UserAuthentication: The return parameters of 'userid', 'height', 'clientID', and 'clientSecret' after verification

The interpretation for the verification:

HS6UserAuthen\_RegisterSuccess: New-user registration succeeded.

HS6UserAuthen\_LoginSuccess: User login succeeded.

HS6UserAuthen\_CombinedSuccess: The user is an iHealth user as well, measurement via SDK has been activated, and the data from the measurement belongs to the user.

HS6UserAuthen\_TrySuccess: Testing without internet connection succeeded.

HS6UserAuthen\_InvalidUserInfo: Userid/clientID/clientSecret verification failed.

HS6UserAuthen\_SDKInvalidateRight: SDK has not been authorized.

HS6UserAuthen\_UserInvalidateRight: User has not been authorized.

HS6UserAuthen\_InternetError: Internet error, verification failed.

The measurement via SDK will be operated in the case of 1-4, and will be terminated if any of 5-8 occurs. The interface needs to be re-called after analyzing the return parameters.

Notice: when a new user registers via SDK, an 'iHealth disclaimer' will pop up automatically, and will require the user to agree in order to continue. SDK applications require an Internet connection; there is 10-day trial period if the SDK cannot connect to the internet, the SDK is fully functional during tryout period, but will be terminated without a working internet connection after 10 days.

binedResult: Action = 1;(1:bind success 2: full user bind failed 3:ts error bind failed 4 other),MAC

Position = 1; (user in weight number)

Status = 1; (1:binding 2: unbind)

TS (time);

binedError: error code

Error code definition:

refer to "error" : HS6 error instruction.

### (3) Unlink HS6

```
-(void)cloudCommandUserBinedQRDeviceWithUser:(HS6User *)tempUser  
deviceID:(NSString *)tempDeviceID  
BlockHS6UserAuthentication:(BlockHS6UserAuthentication)BlockHS6UserAuthenticatio  
n binedResult:(BinedQRDeviceBlock)result  
binedError:(BinedQRDeviceErrorBlock)error;
```

#### **Import parameter:**

tempUser: Properties included: clientID, clientSecret, userID, height.

deviceID: device Mac

#### **Return parameters:**

disBinedResult: Action = 1;(1:bind success 2: full user bind failed 3:ts error bind failed 4 other),MAC

Status = 1; (1:binding 2: unbind)

disBinedError: error code

Error code definition:

refer to "error" : HS6 error instruction.

## 4. Error code instruction

error code:

7: User verify error

101,102,103,104: network error